

Leg Drainage Bags

Continence Aids Product Information Sheet

Hollister Extended Wear Leg Bag – Sterile

(This is not a weekly change leg bag)



Description of Product: Hollister Extended Wear Leg Bag is supplied with adjustable extension tube which attaches directly to a sheath, indwelling catheter or night drainage bag. These bags are not designed to be changed weekly. They are designed to be changed only when the catheter is changed, i.e., every 8 weeks.

NB: Leg straps are supplied with these bags.

Available in:

Product Name	Product Code	Outlet Tap	Tubing Length	Total Capacity	Current MASS Distributor
Hollister Extended Wear Leg Bag	9348	Lever Tap	Adjustable	540mL	Independence
Hollister Extended Wear Leg Bag	9349	Lever Tap	Adjustable	900mL	Independence

Hollister Inview Leg Bag - Sterile



Description of Product: Hollister Inview Leg Drainage Bags attach directly to a sheath, indwelling catheter or night drainage bag. The bag features a non-return valve and fabric backing.

NB: Leg straps are supplied with these bags.

Available in:

Product Name	Product Code	Outlet Tap	Tubing Length	Total Capacity	Current MASS Distributor
Hollister Inview Leg Bag	9621	Lever Tap	10cm	500mL	Independence
Hollister Inview Leg Bag	9624	Lever Tap	50cm	500mL	Independence
Hollister Inview Leg Bag	9631	Lever Tap	10cm	800mL	Independence
Hollister Inview Leg Bag	9632	Lever Tap	50cm	800mL	Independence

Urimaxx (Aaxis) - Sterile



Description of Product: Urimaxx Leg Drainage Bags attach directly to a sheath, indwelling catheter or night drainage bag. The bag features a non-return valve.

NB: Leg straps are supplied with these bags.

Available in:

Product Name	Product Code	Outlet Tap	Tubing Length	Total Capacity	Current MASS Distributor
Urimaxx (Aaxis) Leg bag	10013304	Lever Tap	30cm	500mL	CH2
Urimaxx (Aaxis) Leg bag	10013402	Lever Tap	6cm	750mL	CH2
Urimaxx (Aaxis) Leg bag	10013403	Lever Tap	10cm	750mL	CH2
Urimaxx (Aaxis) Leg bag	10013404	Lever Tap	30cm	750mL	CH2
Urimaxx (Aaxis) Leg bag	10013405	Lever Tap	50cm	750mL	CH2

Conveen Security+ (Coloplast) – Sterile



Description of Product: Conveen Security+ Leg Drainage Bags attach directly to a sheath, indwelling catheter or night drainage bag. The bag features a non-return valve and fabric backing.

NB: Leg straps are supplied with these bags.

Available in:

Product Name	Product Code	Outlet Tap	Tubing Length	Total Capacity	Current MASS Distributor
Conveen Security+ (Coloplast)	5168	Flip-Tap	6cm	500mL	Independence
Conveen Security+ (Coloplast)	51621	Flip-Tap	30cm	500mL	Independence
Conveen Security+ (Coloplast)	51631	Flip-Tap	50cm	500mL	Independence
Conveen Security+ (Coloplast)	210171	Flip-Tap	30cm	750mL	Independence
Conveen Security+ (Coloplast)	210181	Flip-Tap	50cm	750mL	Independence

Bard Uriplan – Sterile



Description of Product: Bard Uriplan Leg Drainage Bags attach directly to a sheath, indwelling catheter or night drainage bag. The bag features a non-return valve and fabric backing.

NB: Leg straps are supplied with these bags.

Available in:

Product Name	Product Code	Outlet Tap	Tubing Length	Total Capacity	Current MASS Distributor
Bard Uriplan	5M	Flip-Tap	10cm	500mL	Independence
Bard Uriplan	5L	Flip-Tap	30cm	500mL	Independence
Bard Uriplan	7M	Flip-Tap	10cm	750mL	Independence
Bard Uriplan	7L	Flip-Tap	30cm	750mL	Independence



© State of Queensland (Queensland Health) 2024
creativecommons.org/licenses/by/3.0/au

For more information contact:

Medical Aids Subsidy Scheme

PO Box 281

Cannon Hill QLD 4170 Australia

tel (07) 3136 3665

MASS-ContinenceAids@health.qld.gov.au

health.qld.gov.au/mass